oclean

The logo features a stylized, thick gray letter 'C' that is open at the top. Below this symbol, the word 'oclean' is written in a lowercase, sans-serif font.

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Thank you for using OcClean Endurance Sonic Electric Toothbrush.

#### Important Safety Tips

If the toothbrush is not used for a long time, the battery will run out due to the power consumption from standby, this can be solved by charging, however, it might take up to 10 minutes of charging in order for the toothbrush to be operational. It is normal for the toothbrush to heat up slightly when it is charging. If the toothbrush becomes extremely hot, immediately stop charging, ensure the toothbrush is turned off and return the product to OcClean's authorized customer service center for inspection. The toothbrush can only be used with the OcClean's brush heads. OcClean will not be able to bear any consequences/responsibilities caused by non-official brush head. The toothbrush is for oral cleaning only. Please do not use it for other purposes. The toothbrush is a personal cleaning item, do not share between more than one user. Please store the toothbrush in a cool and ventilated place and avoid exposure to direct sunlight. Caution for Oral Sensitive Person: The following 3 types of people are not recommended to use this toothbrush unless you have already consulted the professional dentists. 1. The person who has undergone or is under treatment for periodontal surgery. 2. The person with poor fitting dentures, has continuously bleeding gums or any other oral related symptoms. 3. The person whose oral is sensitive. Unless under the supervision of a guardian, this toothbrush is not suitable for people with reduced mental and physical capabilities, those who lack related brushing experience, and children under the age of 12.

Notes: Please do not plug in the charging cable when the toothbrush doesn't need charging.

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### Recommended Gestures for Different Tooth Region

Using Instructions:

please insert the brush head straight onto the toothbrush body. A small gap of about 0.5mm should be left between the head and the main body so the brush head will have space to function properly;

please turn on the toothbrush when it has been put into your mouth with a small amount of soaked toothpaste to prevent toothpaste from splashing;

Hold the toothbrush in a position which is at 45° facing the gum line;

Hold and slowly move the toothbrush in a circular motion, and pass along the gum line. Careful not to press the toothbrush too hard;

Gargle with clean water after brushing;

### Interaction Introduction

**Operation Guide.** Start the toothbrush by short pressing the power button, the toothbrush will automatically begin its countdown of brushing. The toothbrush will pause vibrations after operating for 2 minutes, and shut down automatically after 4 minutes. **Adjustment Period Feature** Oclean Endurance has a function which helps first time electric toothbrush users to adjust its vibration cleaning action. For the first 15 times of operating, the toothbrush enters a state of "adjustment period" and undertake the "soft" intensity level to help the first time user adjusting. Users that feel the vibrating is too weak can skip the adjustment period by long press the power button for 2 seconds and switch to the "Strong" intensity level, the "adjustment period" will not be functioning from that point forward. **Intensity Switch.**

**Wall Mount:** The Wall Mount Accessory contains an adhesive pad at the back, to use, please peel off the protective layer on the back of the cover, and make sure the center decoration is horizontally straight, and positioned the wall mount vertically on the wall. The wall mount can be installed in both upward and downward direction, To remove the wall mount, a blow dryer can be used to heat up the 3M adhesive pad slightly, and pull the wall mount from the installed surface. (carefully not to blow the surface of the wall mount to avoid burnt injuries).

**Charging Process:** When the battery level remaining

is below 5%, its motor will signal by vibrating after brushing section is completed. Please remove the protective cover at the bottom of the toothbrush and insert the Type-C charging cable into the charging port at the bottom of the toothbrush to charge. It will take roughly around 5 hours for the toothbrush to be fully charged. **First Use:** Charging time depends on the battery life demand of user. **Notes:** The battery cannot be replaced by users themselves and can only be trained professionals from Oclean authorized customer service center. Any untrained person should not disassemble this product. If the product was damaged, please contact the Oclean authorized customer service center for maintenance. It is recommended to use Adaptor of 5V/1A for charging, such as the adaptor for iPhone or Tablet.

If the charge current is less than 1A, the charging might take longer.

If the charger voltage is higher than 5V, it could damage the charging circuit.

### FCC WARNING

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. **NOTE:** This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.

Specifications	Product Name	Oclean Endurance Electric Toothbrush	Model Name	E5501	Battery Life	30 Days (*)	Standard	GB39868-2020
	Charging	USB Type-C Charging Cable	Charger	1.8W	Battery Capacity	800 mAh (**)		
	Charging Standard	DC 5V /1A or above	Input	5V / 0.169A	Charging Time	5 Hours (***)		

\* Validated by 2 minute brushing period, twice a day using soft intensity.

\*\* The battery capacity theoretically has a capacity range of  $\pm 50$ mAh.

\*\*\* It is recommended to use a power adapter of 5V/1A

1. Do NOT disassemble, strike, extrude or throw the device into a fire.
2. Do NOT use the toothbrush if its battery is bulging or has expanded.
3. Do NOT place in a high-temperature environment.
4. Do NOT use the device if its battery has been subjected to immersion in water.
5. Before the appliance is discarded, its battery must be removed from it.
6. Before the battery is removed, the appliance must be powered off.
7. The battery shall be disposed of safely.